

# Storing Group Templates

By using the mouse, users can drag *group templates* from the locations listed below into their workplace.

Query Tool Panel (Group label)

The screenshot shows the 'Query Tool' window. At the top, there is a 'Query Name:' field. Below it, there are three columns for group templates. The first column, 'Group 1', is highlighted with a red box and contains a list of medical terms: 'CRP, High Sens. (Cardio) (Group:CPK)', 'CPK (Group:CPK)', 'Troponin-I (Group:TROP-I)', and 'Troponin-T (Group:TROP-T)'. The second column, 'Group 2', is empty. The third column, 'Group 3', is also empty. To the right of the group templates, there is a table titled 'Analysis Types' with the following items:

Analysis Types	
<input type="checkbox"/>	Patient list
<input type="checkbox"/>	Number of pa...
<input type="checkbox"/>	Gender patien...
<input type="checkbox"/>	Vital Status p...
<input type="checkbox"/>	Race patient ...
<input type="checkbox"/>	Age patient b...
<input checked="" type="checkbox"/>	Timeline

At the bottom of the window, there is a 'Run Query Above' button and a 'Patient(s) returned:' field.