

13.4.3.3 Patient in Project Property

Step 6: Check Patient in Project

When a *get* request is sent to the IM cell it has the capability to verify whether or not a patient is part of the project requesting the patient data. By default the property is set to *true*, if you do not wish to check the project then you can set it to *false*.



Note

In addition to setting this property to true, the IM_PROJECT_PATIENTS table has to be populated with the appropriate data in order to check if the patient is in a project. For more information on this table please see the IM Design Document.

1. In the **im.properties** file scroll to the section called **Patient in Project**.
2. Edit the **im.checkPatientInProject** to either *true* or *false*.

EMPI Service

im.checkPatientInProject=**true**

3. Save your changes but do not close the file. You will be continuing to edit the file in the next section.