


Storing Observations

By using the mouse, users can store *observations* such as patient notes from the **Timeline View's Notes Viewer**. Simply place the cursor within the notes viewer and drag to the Workplace view.



The screenshot shows a window titled "Notes Viewer" with a menu bar containing "File" and "Edit". The main content area displays patient information and a table of spirometry results.

Patient Information:

PT#: 00001234 AGE: 32 SEX: M HT: 63.0 in WT: 105.0 lb
PHYSICIAN: abcdef TECH: MAA ASBII

Spirometry Table:

		Pre-Drug*					
		Predicted	Actual	%Pred	Actual	%Pred	%Cha
FVC	(L)	2.58	2.12	82			
FEV1	(L)	1.97	1.51	76			
FEV1/FVC	(%)	77	71	92			
FEF25-75%	(L/S)	1.79	1.09	61			
FEFmax	(L/S)	5.16	2.70	52			
TET	(SEC)		9.66				