

## 3.7.1 Working Directory

### Step1: Change Working Directory for Metadata Installation

During the **Metadata Installation** we will be working with the **Metadata directory**. Within this directory are the files we will be editing to define the database properties, create the IM tables and load the sample data provided with the data package.

To change your working directory to the Metadata directory enter the following command:

#### Change Directory Command

```
CD YOUR_I2B2_DATA_DIRedu.harvard.i2b2.data\Release_1-7\NewInstall\Metadata
```



#### Important

The change directory command shown above is written for Windows. The command itself is the same for both Linux and Windows with the only exceptions being the slashes in the path name and Linux does not include the "drive".

- *Linux Path:* use a forward slash.
- *Windows Path:* use a backward slash.

#### EXAMPLES:

**Windows Command:** CD C:\opt\data\edu.harvard.i2b2.data\Release\_1-7\NewInstall\Metadata

**Linux Command:** CD /opt/data/edu.harvard.i2b2.data/Release\_1-7/NewInstall/Metadata