## Quick Upgrade Guide for i2b2 Web Client 1.8.x

- 1. Download the new v1.8.x Web Client release.
- 2. Create a new hosting directory (http://hostname/legacy) in Apache that will contain your currently existing deployment.
- 3. Move the files for the existing web client deployment into this new hosting directory
- 4. Copy the new v1.8 web client into the root deployment directory (http://hostname/).
- 5. Edit proxy.php and setup the whitelist and blacklist entries.
- 6. Edit i2b2\_config\_domains.js file and change the "urlProxy" entry to "/proxy.php".
  - a. Copy the i2b2 domains from your legacy deployment into the i2b2\_config\_domains.js file. See "Apache Install Guide for i2b2 Release v1. 8+" section "I2B2 CONFIGURATION".
  - b. Configure any SAML (Shibboleth) authentication within the i2b2 domains if needed.
- 7. Download any additional "modern" plugins and move into the /plugins directory.
- Edit /plugins/plugins.json file and add the new plugins. See "Apache Install Guide for i2b2 Release v1.8+" section "CONFIGURATION OF MODERN PLUGINS".
- 9. Copy into /js-i2b2/cells/LEGACYPLUGIN/legacy\_plugin directory any legacy plugins from the /legacy deployment.
- 10. Register the copied plugins (from the previous step) into the following file:
  - /js-i2b2/cells/LEGACYPLUGIN/legacy\_plugin/js-i2b2/i2b2\_loader.js
  - See "Install Guide for i2b2 Release v1.8+" section "CONFIGURATION OF LEGACY PLUGINS" for examples and more information.
- 11. Edit the login screen template:
  - a. Set the URL address for the "classic i2b2" link to "/legacy".
  - b. Set the URL address for the "Forgot password" link if you have a custom page to reset a user's password. Remove this linkfrom the HTML template if you do not offer this service.
  - c. Add any additional SAML authentication service links if you use something besides Shibboleth or Okta.
- 12. Apply any logos, color schemes or additional branding to the template files. See the document titled "Applying Branding Styles to the i2b2 Release v1.8+" for more details.